STERYLAB PRODUCT RANGE

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**THE BONE MARROW TREPINE BIOPSY NEEDLE WHICH DOES
NOT NECESSITATE LUXATION OR ADDITIONAL MOVEMENTS**


best®
L I S A S

**LESS INVASIVE SPECIMEN
ACQUIREMENT SYSTEM**

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HOW THE NOLOSE® SYSTEM WORKS

The cannulas TREPAN® cutting profile allows the attainment of uniform intact osteo-medullary cylinders. The inner sample retaining system requires just a moderate luxatory movement.

In case of soft or critical bones, or everytime luxatory manoeuvre is to be avoided, the use combined with the SAFE-LOCK® accessory (included), will guarantee success in specimen acquirement.

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9G (3,5 mm) / 11G (3,0 mm)

8G (4,0 mm) / 11G (3,0 mm)

Ø Frustolo biottico / Ø of specimen core

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THE BONE MARROW TREPHINE BIOPSY NEEDLE WHICH DOES NOT NECESSITATE LUXATION OR ADDITIONAL MOVEMENTS

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 The advanced design of BEST-LISAS® will guarantee the following advantages:

► **Ease of use and guaranteed specimen collection**

The unique NOLOSE® core retention system is easy to use even for the inexperienced for the most common cases. For soft bones or critical pathologies the use of BEST-LISAS® needle combined with the accessory SAFE-LOCK® will always guarantee success in the specimen acquisition. NO MORE LOST CORES!

► **A considerable reduction in pain**

The smaller Gauge of BEST-LISAS® will give a reduction in puncture area over a traditional Jamshidi needle (see comparative chart) with the same specimen volume and, combined with no need to luxate, renders the procedure particularly pain free.

► **A more compact and abundant specimen**

The TREPAN® crown edge with its advanced geometry will guarantee an intact and larger specimen core even in the softest bone with less pain to the patient and ease of use for the physician.

► **A simpler and more rapid procedure**

BEST-LISAS® requires no luxatory movement to retain the core specimen.

► Colour coded for easier needle size (gauge) determination

► Sterile, single use: reduction of infection hazard.

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STERYLAB S.r.l. • VIA MAGENTA 77/6 • 20017 RHO (MILANO) • ITALY